


# $\mu$ Contact Printers

Instruments for Microcontact Printing and NIL




Create 2D and 3D Structures on the Micro- and Nanoscale





## μCP Platforms



μCP4.1 with all items including fluidic/pneumatic control box and external UV lamp (OmniCure)



Two-level μCP6.1 with fast linear XY-motors. The lower Y-axis contains a UV collimator for precise and homogeneous irradiation especially for UV-NIL (see back page), the upper one the stamping unit. The flexible design allows easy reconfiguration.

### The making of micro- and nanostructures in 2D...

Few affordable techniques exist to produce micro-patterns on a surface. Microarray spotting (e.g. with the GeSiM Nano-Plotter) generates resolutions around 100 μm. **Microcontact printing (μCP)**, on the other hand, can work on the nanoscale: a stamp of an elastomeric polymer like PDMS (silicone), cast on a micro- or nanostructured master, transfers molecules to an even surface. Chemicals, self-assembling monolayers (SAMs), biomolecules, nanoparticles, beads and cells can be printed with unlimited complexity, which is why μCP is getting increasingly popular in the life sciences, e.g. to study the influence of matrix proteins or growth factors on growth, differentiation and movement of cells.

### ...and 3D

**Nanoimprint lithography (NIL)** has become a hot topic in microelectronics and cell biology. Here the 3D stamp structure is replicated in a polymer at elevated temperature and fixed by cooling or UV cross-linking.



Both μCP and NIL can be performed on the automatic GeSiM μContactPrinters or the semi-automatic μCP Core, making both methods completely reproducible. The printed structures can be tiled together to a larger one via step & repeat, especially when doing UV-NIL with a silicon mask in the stamp to exactly define the area that is cured ('double-side patterning'). For higher throughput and more exactly defined UV curing using a collimator on the lower level, the μCP6.1 exists. The systems speak G-code and their repeating accuracy in X/Y is better than 5 μm.

### Stamps

Shipment includes a casting station that lets you manufacture your own stamps at reasonable costs. Customized, teflonised silicon masters with structure sizes as small as 100 nm can be ordered from GeSiM. Stamps are usually made of PDMS, except for NIL structures < 500 nm with high aspect ratio that require stiffer materials such as PFPE.

Stamps are picked up from a rack, allowing to combine different stamps during a run, e.g. when inked with different proteins. The accurate pneumatic actuation of the stamp (for NIL: high overpressure while stamp is down and de-tooling using vacuum pulses while lifting up; GeSiM specialty!) is fundamental for reproducible prints across the entire area.

An integrated microscope and the rotating print head help to accurately align stamp and substrate (picture), either manually or via image processing.



μCP Core, here with large stamp and UV optical fibre, on an inverted microscope (control box, temperature control, UV lamp and tubes for compressed air/vacuum not shown)


### Platforms

All GeSiM μContact Printers feature a proprietary stamping technology based on an accurate Z-drive and pneumatic actuation during printing. It is precise and reproducible and does not need expensive granite beds and vibration damping tables.


- **μCP4.x:** The smallest platform with XYZ stage, with camera, stamp holder and UV lens head on a single level. Max. stamp size 20 mm × 20 mm.
- **μCP6.x:** A μCP with linear motors. A UV source with collimator and photo mask on the lower level irradiates the transparent substrate on the upper level. The alignment of the exposition frame with the stamping area allows 'chessboard'-like multiple imprints on one target (up to 6" = 152 mm large) with zero spacing in between.
- **μCP6.x/E:** Like μCP6.x, with twice the work area
- **μCP Core:** Semi-automatic, budget μCP device consisting of μCP Z-drive and adjustable substrate holder. It provides the same stamping technology as the other instruments, but is placed on an inverted microscope (not included).

### Patterns for cell growth

One can modulate cell growth and hence differentiation by  $\mu$ -contact printing of patches of extracellular matrix proteins. Here, 100  $\mu$ m large squares of fibronectin were microcontact-printed with a PDMS stamp onto a cell-repellent, plasma-polymerized polyethylene glycol (PEG), and human umbilical cord blood neural stem cells (HUCB-NSCs) were grown on these squares to study the influence of cell motility on differentiation, as cells adhere only to the fibronectin squares. After culturing, cells were fixed and immunostained. Courtesy of EC Joint Research Centre, Nanobiosciences Unit, Ispra, Italy (Dora Mehn)



HUCB-NSC cells on 100  $\mu$ m fibronectin squares, fluorescence image. Green: beta-tubulin-III neuronal marker; red: GFAP; blue: cell nuclei.



Master structure for stamp casting, made by 2-photon lithography using a stiff material.

Micro-needle array from Norland NOA63 optical adhesive, printed using a PDMS stamp.

### Micro-needle arrays

Micro-needle arrays are gaining interest for non-invasive, convenient transdermal drug delivery. Micro-needles consist of bioresorbable polymers blended with deliverables or of solid material (with reservoirs or coated with active ingredients). NIL allows better spatial resolution than 3D-SLA, at low operational costs and sizes  $\leq 1$  mm, thus allowing quick processes using a wide range of materials. GeSiM's NIL technology enables micro-needles with grooves, indentations and channels. Subsequent drug deposition onto them can be done with GeSiM piezo pipettes, either directly in the  $\mu$ CP or on a GeSiM Nano-Plotter.


High aspect ratios (here: height 1.2 mm, radius 0.2 mm) are achieved by evacuating the space around the PDMS stamp while it is in contact with the surface.

### Nanoplasmonic sensors


Surface plasmon resonance (SPR) sensors enable the label-free real-time quantification of biomolecules. Classical SPR utilizes a gold film on a prism. To scale this down, a nanostructured gold layer was used to excite plasmons without bulky optics. The resulting sensor was used to measure pollutants in wastewater. Central to this method was the reproducible fabrication of nanostructures by UV-NIL using a GeSiM  $\mu$ CP 4.1 and a UV-cross-linkable photo-resist, followed by gold coating. Ref.: Steinke, N. et al., *Sens. Act. B* 254, 2018, 749-754

Similarly, ethanol formation during fermentation can be measured with an ethanol-sensitive hydrogel with nanoplasmonic surface.

Ref.: Kroh, C. et al., *J. Sens. Sens. Syst.* 7, 2018, 51-55



Nano-patterns repeatedly printed on a glass wafer. Small picture: REM image of the nano-pillars; the pitch is 450 nm.



Picowells, 25  $\mu$ m diameter (large picture), and spheroids in 250  $\mu$ m wide honeycomb wells, both formed via NIL

### Picolitre wells for cell biology

In cooperation with M. Deutsch's lab from Bar-Ilan University and others, pico-well arrays were printed by UV-NIL in UV-curable adhesives such as NOA81 to measuring the activity of individual cells e.g. via multi-parametric enzymatic assays. One application was toxicology testing of nanoparticles that the GeSiM Nano-Plotter had spotted into the picowells, another the formation and study of organoids from cancer cells in honeycomb structures, with and without electrodes. Courtesy of Motti Deutsch, Bar-Ilan University.


Ref. (selection): Afrimson, E. et al., *BMC Cell Biology* 11, 2010, 83

Markovitz-Bishitz, Y. et al., *Biomaterials* 31, 2010, 8436-8444

Zurgil, N. et al., *Lab Chip* 14, 2014, 2226-2239


### In vitro neurotoxicity assay development

To direct neuronal growth to certain sites, 100  $\mu$ m wells above microelectrodes and channels between the wells were formed by UV-NIL. The stamp was aligned with the electrodes and a PEG-based hydrogel polymerized by 60 s UV exposition (left picture). The right picture shows rat neurons grown in such wells/channels (green:  $\beta$ -tubulin III neuronal marker, white bar 50  $\mu$ m). Courtesy of EC Joint Research Centre, Nanobiosciences Unit, Ispra, Italy (Dora Mehn, Jakub Nowak).






### The revolutionary stamping technology of the GeSiM $\mu$ ContactPrinters



GeSiM  $\mu$ CP stamping principle: the stamp membrane is bulged out while contacting the surface, resulting in defined pressure and even transfer. Demoulding after NIL utilizes controlled pulses of vacuum while moving up; this ensures that structures with high aspect ratios don't get destroyed.

### Stamp casting




Casting stations for small and large stamps. Max. stamp size is 20x20 (mmxmm).

2022


### 'Double side patterning' / 'chessboard-like' UV-NIL

The cured region in UV-NIL is normally larger than the stamp region. Placing a silicon mask in the stamp, however, restricts the cross-linked (cured) area to the mask area; this allows to exactly place NIL structures side by side in a chequerboard-like manner:



Double-side patterning with the  $\mu$ CP4.x. The silicon mask limits the size of the cross-linked area under the stamp.


In the  $\mu$ CP6.x, a **collimator** on the lower level produces parallel UV light beams that shine through a photo mask before they reach the transparent substrate on the upper level. Not only is the cross-linked area more exactly defined, you can also insert a *patternd* mask to use two structuring methods at once, namely photolithography and stamping. This can generate quite complex designs, e.g. a nanostructure in a microfluidic channel, in a single step.



UV collimator (cover removed) on the lower Y-axis of a  $\mu$ CP6.1. It allows to mount a photo mask and very precise adjustment via micrometre screws.

### Tools on the work deck and on the multi-Z tool head

The stamping unit is assisted by other tools. A typical configuration of the  $\mu$ CP4.x is shown here:



Work deck tools: stamp rack, spin coater (e.g. for semi-dry inking, with stamp height sensor), heatable/coolable (up to 100 °C) holder for 4-inch substrates. In the rear, hidden: wash/dry stations, tip aligner, stroboscope. Tools on the print head: stamping unit, piezo pipette and substrate height sensor (on same Z-axis), heatable cartridge holder for pastes/glues, especially for NIL.

Select from a plethora of accessories, also from other GeSiM Robotics systems, i.e. BioScaffolder (3D bioprinter) or the BSYS (pipetting machine for chemical synthesis), so please see their brochures for details:

- Integrated dust-protecting lid, making you independent of biological safety cabinets
- Work plate: coolable/heatable microtitre plate holders, 'double-side patterning' stamps (see left, with casting station), holders for slides, glass vials and microcentrifuge tubes
- Print head: various dispensers, e.g. cartridge cooler (5...100 °C), piezo-driven powder dispensers for  $\mu$ g quantities, pipettors with disposable plastic tips, pH titration tip, plasma pen etc.

Gesellschaft für Silizium-Mikrosysteme mbH  
Bautzner Landstraße 45  
01454 Radeberg, Germany  
Tel. +49-351-2695 322  
Fax +49-351-2695 320  
contact@gesim.de  
www.gesim.de



For more information or to see the latest product changes and enhancements, please visit our website or call.

